

## OmpC Protein, Klebsiella pneumoniae, Recombinant (His &amp; Myc)

## General Information

Synonyms:	Outer membrane porin C; Outer membrane protein C; Porin OmpC; ompC; Porin ompk36; ompK36
Protein Construction:	22-363 aa
Species:	Klebsiella pneumoniae
Expression Host:	E. coli
Accession:	Q48473
Molecular Weight:	45.0 kDa (predicted)
AA Sequence:	AEIYNKDGKLDLYGKIDGLHYFSDDKDVDGDQTYMRLGVKGETQINDQLTGYGQWEYNVQANNTSSSD QAWTRLAFAGLKFGDAGSFDYGRNYGVVYDVTSWTDVLPFEGGDTYGSDFLQSRANGVATYRNSDFFGL VDGLNFALQYQKNGSVSGEGATNNGRGALKQNGDGFSTVYDIFDGISAGFAYANSKRTDDQNQLLGE GDHAETYTGGLKYDANNIYLATQYTQTYNATRAGSLGFANKAQNFVAAQYQDFGLRPSVAYLQSKGKDL NGYGDQDILKYVDVGATYYFNKNMSTYVDYKINLLDDNSFTRSAGISTDDVVALGLVYQF

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100  $\mu$ g/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

## Protein Background

## A DRUG SCREENING EXPERT

---

Forms pores that allow passive diffusion of small molecules across the outer membrane. In *K.pneumoniae* it has been shown to bind the C1Q component and activate the classical pathway of the complement system.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481