

Anti-Myoglobin Antibody (2Q412)

Product Details

Ig Type:	IgG
Reactivity:	Human,Rat,Mouse
Molecular Weight:	Theoretical: 17 kDa. Actual: 17 kDa.
Clone:	2Q412
Purification:	Protein G purified

Applications

1. Paraformaldehyde-fixed, paraffin embedded (Rat skeletal muscle); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (Myoglobin) Polyclonal Antibody, Unconjugated (TMAB-09174) at 1: 500 overnight at 4°C, followed by a conjugated secondary for 20 minutes and DAB staining.
2. Paraformaldehyde-fixed, paraffin embedded (human skeletal muscle); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (Myoglobin) Monoclonal Antibody, Unconjugated (TMAB-09174) at 1: 200 overnight at 4°C, followed by operating according to SP Kit (Mouse) instructions and DAB staining.
3. Paraformaldehyde-fixed, paraffin embedded (mouse skeletal muscle); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30 min; Antibody incubation with (Myoglobin) Monoclonal Antibody, Unconjugated (TMAB-09174) at 1: 200 overnight at 4°C, followed by operating according to SP Kit (Mouse) instructions and DAB staining.
4. Sample:
Lane 1: Rat Muscle tissue lysates
Primary: Anti-Myoglobin (TMAB-09174) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution
Predicted band size: 17 kDa
Observed band size: 17 kDa
5. (Negative control) paraformaldehyde-fixed, paraffin embedded Rat Liver; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Myoglobin Monoclonal Antibody, Unconjugated (TMAB-09174) at 1: 200 overnight at 4°C, followed by conjugation to the AffiniPure Goat Anti-Mouse IgG H&L-HRP and DAB staining.
6. (Negative control) paraformaldehyde-fixed, paraffin embedded Human Liver; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Myoglobin Monoclonal Antibody, Unconjugated (TMAB-09174) at 1: 200 overnight at 4°C, followed by conjugation to the AffiniPure Goat Anti-Mouse IgG H&L-HRP and DAB staining.
7. Paraformaldehyde-fixed, paraffin embedded Rat Heart; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Myoglobin Monoclonal Antibody, Unconjugated (TMAB-09174) at 1: 200 overnight at 4°C, followed by conjugation to the AffiniPure Goat Anti-Mouse IgG H&L-HRP and DAB staining.
8. Paraformaldehyde-fixed, paraffin embedded Human Heart; Antigen retrieval by boiling in

Verified Activity:

A DRUG SCREENING EXPERT

sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Myoglobin Monoclonal Antibody, Unconjugated (TMAB-09174) at 1: 200 overnight at 4°C, followed by conjugation to the AffiniPure Goat Anti-Mouse IgG H&L-HRP and DAB staining.

Application: WB,IHC-P,IHC-Fr,IF

Recommended WB: 1:500-2000; IHC-P: 1:100-500; IHC-Fr: 1:100-500; IF: 1:100-500

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Purified Human Myoglobin (Natural Human Myoglobin protein)

Antigen Species: Human

Gene ID: 4151

Uniprot ID: P02144

Research Background

Myoglobin is a small heme containing protein responsible for the oxygen deposition in muscle tissues. Only one form of myoglobin is expressed in cardiac and skeletal muscles. Myoglobin is known as a marker of myocardial damage and it has been used for more than three decades. Nowadays it still is very commonly used in clinical practice as an early marker of AMI. It appears in patient's blood 1 to 3 hours after onset of the symptoms, reaching peak level within 8 to 12 hours. Myoglobin is not so cardiac specific as cTnI or cTnT. Because of high myoglobin concentration in skeletal muscle tissue, even minor skeletal muscle injury results in the significant increase of myoglobin concentration in blood. Thus myoglobin is used together with cTnI or cTnT in clinical practise for better specificity in AMI diagnosis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481