Data Sheet (Cat.No.T71120)



AR-13503

Chemical Properties


CAS No.: 1254032-16-0

Formula: C19H19N3O2

Molecular Weight: 321.37

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description AR-13503, is an metabolite of Netarsudil (AR-11324), is a Rho-associated protein kinase

inhibitor. Netarsudil is potential useful for treating glaucoma and/or reducing

intraocular pressure.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1117 mL	15.5584 mL	31.1168 mL
5 mM	0.6223 mL	3.1117 mL	6.2234 mL
10 mM	0.3112 mL	1.5558 mL	3.1117 mL
50 mM	0.0622 mL	0.3112 mL	0.6223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com