


(6-methoxyppyridin-3-yl)methanol

Chemical Properties

CAS No. :	58584-63-7
Formula:	C7H9NO2
Molecular Weight:	139.154
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(6-methoxyppyridin-3-yl)methanol, with CAS No. 58584-63-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (6-methoxyppyridin-3-yl)methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Solubility Information

Solubility	DMSO: 27.83 mg/mL (199.99 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.1865 mL	35.9324 mL	71.8649 mL
5 mM	1.4373 mL	7.1865 mL	14.373 mL
10 mM	0.7186 mL	3.5932 mL	7.1865 mL
50 mM	0.1437 mL	0.7186 mL	1.4373 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481