


## Thonningianin B

## Chemical Properties

CAS No. :	271579-12-5
Formula:	C35H30O17
Molecular Weight:	722.608
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Thonningianin B is a natural product from Penthorum chinense Pursh
Targets(IC50)	Others, Autophagy

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3839 mL	6.9194 mL	13.8387 mL
5 mM	0.2768 mL	1.3839 mL	2.7677 mL
10 mM	0.1384 mL	0.6919 mL	1.3839 mL
50 mM	0.0277 mL	0.1384 mL	0.2768 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481