


## Iptacopan hydrochloride hydrate

## Chemical Properties

CAS No. :	2447007-60-3
Formula:	C <sub>25</sub> H <sub>33</sub> ClN <sub>2</sub> O <sub>5</sub>
Molecular Weight:	476.99
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Iptacopan hydrochloride hydrate is a highly potent and highly selective factor B inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0965 mL	10.4824 mL	20.9648 mL
5 mM	0.4193 mL	2.0965 mL	4.193 mL
10 mM	0.2096 mL	1.0482 mL	2.0965 mL
50 mM	0.0419 mL	0.2096 mL	0.4193 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481