


## IMP2-IN-2

## Chemical Properties

CAS No. :	1448353-39-6
Formula:	C <sub>22</sub> H <sub>19</sub> F <sub>3</sub> N <sub>2</sub> O <sub>3</sub> S
Molecular Weight:	448.46
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	IMP2-IN-2 is a potent and selective inhibitor of IMP2 with IC <sub>50</sub> values of 120.9 μM and 236.7 μM for the inhibition of IMP2 interaction with RNA_A and RNA_B, respectively. IMP2-IN-2 can be used in cancer research.
Targets(IC <sub>50</sub> )	Others,IGF-1R

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2299 mL	11.1493 mL	22.2985 mL
5 mM	0.446 mL	2.2299 mL	4.4597 mL
10 mM	0.223 mL	1.1149 mL	2.2299 mL
50 mM	0.0446 mL	0.223 mL	0.446 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481