


2-Chloroaniline hydrochloride

Chemical Properties

CAS No. :	137-04-2
Formula:	C ₆ H ₆ ClN.ClH
Molecular Weight:	164.03
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

HCl



Biological Description

Description	2-Chloroaniline hydrochloride, with CAS No. 137-04-2, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-Chloroaniline hydrochloride provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0964 mL	30.4822 mL	60.9645 mL
5 mM	1.2193 mL	6.0964 mL	12.1929 mL
10 mM	0.6096 mL	3.0482 mL	6.0964 mL
50 mM	0.1219 mL	0.6096 mL	1.2193 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481