


Mol4

Chemical Properties

CAS No. :	512800-53-2
Formula:	C30H20ClNO5S
Molecular Weight:	542
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mol4 is a novel BCL-2 inhibitor (IC50: 153.3 µM).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.845 mL	9.2251 mL	18.4502 mL
5 mM	0.369 mL	1.845 mL	3.690 mL
10 mM	0.1845 mL	0.9225 mL	1.845 mL
50 mM	0.0369 mL	0.1845 mL	0.369 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481