


## Famiraprinium chloride

## Chemical Properties

CAS No. : 96440-63-0  
 Formula: C<sub>15</sub>H<sub>18</sub>ClN<sub>3</sub>O<sub>2</sub>  
 Molecular Weight: 307.78  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Famiraprinium chloride is a competitive antagonist at the GABA receptor site.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2491 mL	16.2454 mL	32.4907 mL
5 mM	0.6498 mL	3.2491 mL	6.4981 mL
10 mM	0.3249 mL	1.6245 mL	3.2491 mL
50 mM	0.065 mL	0.3249 mL	0.6498 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481