


Vinrosidine sulfate

Chemical Properties

CAS No. :	18556-44-0
Formula:	C ₄₆ H ₆₀ N ₄ O ₁₃ S
Molecular Weight:	909.06
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Vinrosidine sulfate is an antineoplastic.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.100 mL	5.5002 mL	11.0004 mL
5 mM	0.220 mL	1.100 mL	2.2001 mL
10 mM	0.110 mL	0.550 mL	1.100 mL
50 mM	0.022 mL	0.110 mL	0.220 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Fyfe MJ, Lofffield S, Goldman ID. A reduction in energy-dependent amino acid transport by microtubular inhibitors in Ehrlich ascites tumor cells. *J Cell Physiol.* 1975 Oct;86(2 Pt 1):201-11. PubMed PMID: 1194361.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481