

Compound T128744

Chemical Properties

CAS No. :

Formula: C35H62O15

Molecular Weight: 722.87

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3-methyl-1-[(3-methyl-1-[(3-methyl-1-((3-methyl-1-oxo-1-[(2,3,4,5-tetrahydroxy)pentan-2-yl]oxy)-1-oxopentan-2-yl]oxy)-1-oxopentan-2-yl]oxy)-1-oxopentan-2-yl]oxy)-1-oxopentan-2-yl 2-hydroxy-3-methylpentanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3834 mL	6.9169 mL	13.8337 mL
5 mM	0.2767 mL	1.3834 mL	2.7667 mL
10 mM	0.1383 mL	0.6917 mL	1.3834 mL
50 mM	0.0277 mL	0.1383 mL	0.2767 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481