

ART615

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>21</sub>F<sub>3</sub>N<sub>4</sub>O<sub>2</sub>

Molecular Weight: 418.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ART615 is a related isomer of ART558. ART615 inhibits Polθ by <10% at 12 μM and is able to act as a control for ART558 (IC <sub>50</sub> :7.9 nM).
Targets(IC <sub>50</sub> )	Others,DNA/RNA Synthesis

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.390 mL	11.950 mL	23.900 mL
5 mM	0.478 mL	2.390 mL	4.780 mL
10 mM	0.239 mL	1.195 mL	2.390 mL
50 mM	0.0478 mL	0.239 mL	0.478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481