

2-[2-(3,4-dihydroxybenzoyloxy)-4,6-dihydroxyphenyl]acetic acid

Chemical Properties


CAS No. :

Formula: C₁₅H₁₂O₈

Molecular Weight: 320.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-[2-(3,4-dihydroxybenzoyloxy)-4,6-dihydroxyphenyl]acetic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1226 mL	15.6128 mL	31.2256 mL
5 mM	0.6245 mL	3.1226 mL	6.2451 mL
10 mM	0.3123 mL	1.5613 mL	3.1226 mL
50 mM	0.0625 mL	0.3123 mL	0.6245 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481