


S 332

Chemical Properties

CAS No. : 53297-82-8
 Formula: C₁₆H₁₇N₃O₃S
 Molecular Weight: 331.39
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	S 332 is an adrenergic beta3-agonist that is more potent in increasing the urinary storage volume than clenbuterol and less potent in inhibiting the contractile force of urinary bladder at micturition reflex than clenbuterol.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0176 mL	15.088 mL	30.1759 mL
5 mM	0.6035 mL	3.0176 mL	6.0352 mL
10 mM	0.3018 mL	1.5088 mL	3.0176 mL
50 mM	0.0604 mL	0.3018 mL	0.6035 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481