

CGS 24891

Chemical Properties


CAS No. : 134823-10-2

Formula: C₁₈H₁₇FN₂O₄

Molecular Weight: 344.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CGS 24891 is a 5-lipoxygenase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9041 mL	14.5205 mL	29.0411 mL
5 mM	0.5808 mL	2.9041 mL	5.8082 mL
10 mM	0.2904 mL	1.4521 mL	2.9041 mL
50 mM	0.0581 mL	0.2904 mL	0.5808 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481