

Compound N016-0020

Chemical Properties


CAS No. :

Formula: C70H105N3O24S3

Molecular Weight: 1468.79

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N016-0020 is a useful organic compound for research related to life sciences and the catalog number is T131348.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6808 mL	3.4042 mL	6.8083 mL
5 mM	0.1362 mL	0.6808 mL	1.3617 mL
10 mM	0.0681 mL	0.3404 mL	0.6808 mL
50 mM	0.0136 mL	0.0681 mL	0.1362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481