


Compound FL0161

Chemical Properties

CAS No. :	95533-25-8
Formula:	C ₉ H ₁₇ NO ₂
Molecular Weight:	171.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound FL0161, with CAS No. 95533-25-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Compound FL0161 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8398 mL	29.1988 mL	58.3976 mL
5 mM	1.168 mL	5.8398 mL	11.6795 mL
10 mM	0.584 mL	2.9199 mL	5.8398 mL
50 mM	0.1168 mL	0.584 mL	1.168 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481