

Compound N037-0036

Chemical Properties


CAS No. :

Formula: C₂₇H₄₅NO₃

Molecular Weight: 431.661

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N037-0036 is a useful organic compound for research related to life sciences and the catalog number is T131287.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3166 mL	11.5832 mL	23.1664 mL
5 mM	0.4633 mL	2.3166 mL	4.6333 mL
10 mM	0.2317 mL	1.1583 mL	2.3166 mL
50 mM	0.0463 mL	0.2317 mL	0.4633 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481