


RO-11-7330

## Chemical Properties

CAS No. : 64301-64-0  
 Formula: C<sub>20</sub>H<sub>17</sub>ClF<sub>3</sub>NS  
 Molecular Weight: 395.87  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	RO-11-7330 is a bioactive chemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5261 mL	12.6304 mL	25.2608 mL
5 mM	0.5052 mL	2.5261 mL	5.0522 mL
10 mM	0.2526 mL	1.263 mL	2.5261 mL
50 mM	0.0505 mL	0.2526 mL	0.5052 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

RO-11-7330 RO-11-7330 Free Base

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481