

## Compound T128219

## Chemical Properties


CAS No. :

Formula: C63H100O30

Molecular Weight: 1337.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-(((2S,3S,4S,5R,6S)-6-(((2S,3R,4S,5S)-2-(((4aS,6aS,6bR,9R,10R,11S,12aR,14bS)-10-(((2R,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)oxan-2-yl]oxy)-11-hydroxy-9-(hydroxymethyl)-2,2,6a,6b,9,12a-hexamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carbonyloxy)-4,5-dihydroxyoxan-3-yl]oxy)-5-hydroxy-2-methyl-4-(((2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy)oxan-3-yl]oxy)-3-hydroxy-3-methyl-5-oxopentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7477 mL	3.7384 mL	7.4769 mL
5 mM	0.1495 mL	0.7477 mL	1.4954 mL
10 mM	0.0748 mL	0.3738 mL	0.7477 mL
50 mM	0.015 mL	0.0748 mL	0.1495 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481