

U 72745G

Chemical Properties


CAS No. : 378-61-0

Formula: C₂₂H₂₉F₀O₄

Molecular Weight: 376.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	U 72745G is an angiogenesis factor antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6563 mL	13.2816 mL	26.5632 mL
5 mM	0.5313 mL	2.6563 mL	5.3126 mL
10 mM	0.2656 mL	1.3282 mL	2.6563 mL
50 mM	0.0531 mL	0.2656 mL	0.5313 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481