

Compound CSC000682698

## Chemical Properties


CAS No. : 123530-67-6

Formula: C<sub>10</sub>H<sub>7</sub>NO<sub>2</sub>S

Molecular Weight: 205.233

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound CSC000682698 is a high-purity biochemical reagent with CAS number 123530-67-6, suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8726 mL	24.3629 mL	48.7258 mL
5 mM	0.9745 mL	4.8726 mL	9.7452 mL
10 mM	0.4873 mL	2.4363 mL	4.8726 mL
50 mM	0.0975 mL	0.4873 mL	0.9745 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481