


Serazapine

Chemical Properties

CAS No. :	115313-22-9
Formula:	C ₂₂ H ₂₃ N ₃ O ₂
Molecular Weight:	361.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Serazapine is a highly specific serotonin (5-HT ₂) binding inhibitor with anxiolytic activity for the treatment of anxiety disorders.
Targets(IC ₅₀)	5-HT Receptor

Solubility Information

Solubility	DMSO: 50 mg/mL (138.34 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7667 mL	13.8336 mL	27.6671 mL
5 mM	0.5533 mL	2.7667 mL	5.5334 mL
10 mM	0.2767 mL	1.3834 mL	2.7667 mL
50 mM	0.0553 mL	0.2767 mL	0.5533 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Katz RJ, et al. Serotonergic (5-HT₂) mediation of anxiety-therapeutic effects of serazapine in generalized anxiety disorder. Biol Psychiatry. 1993;34(1-2):41-44.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481