

DB03417

## Chemical Properties


CAS No. : 1027198-20-4

Formula: C<sub>18</sub>H<sub>33</sub>N<sub>7</sub>O<sub>5</sub>

Molecular Weight: 427.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	DB03417 is a compound listed in web page: <a href="https://go.drugbank.com/drugs/DB03417">https://go.drugbank.com/drugs/DB03417</a> . DB03417 is a potential tryptase inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3392 mL	11.6959 mL	23.3918 mL
5 mM	0.4678 mL	2.3392 mL	4.6784 mL
10 mM	0.2339 mL	1.1696 mL	2.3392 mL
50 mM	0.0468 mL	0.2339 mL	0.4678 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481