

Ethyl tetrahydropyran-4-yl-acetate

Chemical Properties

CAS No. :	103260-44-2
Formula:	C ₉ H ₁₆ O ₃
Molecular Weight:	172.221
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Ethyl tetrahydropyran-4-yl-acetate ,with CAS No. 103260-44-2, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Ethyl tetrahydropyran-4-yl-acetate provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8065 mL	29.0326 mL	58.0653 mL
5 mM	1.1613 mL	5.8065 mL	11.6131 mL
10 mM	0.5807 mL	2.9033 mL	5.8065 mL
50 mM	0.1161 mL	0.5807 mL	1.1613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481