

PPAR α agonist 1

Chemical Properties

CAS No. :

Formula: C₂₇H₃₄O₄

Molecular Weight: 422.56

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PPAR α agonist 1 is a complete and potent PPAR α agonist.
Targets(IC50)	Others,PPAR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3665 mL	11.8326 mL	23.6653 mL
5 mM	0.4733 mL	2.3665 mL	4.7331 mL
10 mM	0.2367 mL	1.1833 mL	2.3665 mL
50 mM	0.0473 mL	0.2367 mL	0.4733 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481