

Compound N094-0027

Chemical Properties


CAS No. :

Formula: C₂₄H₃₉NO₅

Molecular Weight: 421.578

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N094-0027 is a useful organic compound for research related to life sciences and the catalog number is T131157.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.372 mL	11.8601 mL	23.7203 mL
5 mM	0.4744 mL	2.372 mL	4.7441 mL
10 mM	0.2372 mL	1.186 mL	2.372 mL
50 mM	0.0474 mL	0.2372 mL	0.4744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481