


(S)-1-Benzylpyrrolidin-3-ol

Chemical Properties

CAS No. :	101385-90-4
Formula:	C11H15NO
Molecular Weight:	177.247
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-1-Benzylpyrrolidin-3-ol is a small molecule compound that can be used to synthesize antihypertensive active molecules.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6417 mL	28.2087 mL	56.4175 mL
5 mM	1.1283 mL	5.6417 mL	11.2835 mL
10 mM	0.5642 mL	2.8209 mL	5.6417 mL
50 mM	0.1128 mL	0.5642 mL	1.1283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481