

Podoblastin A (D-mannitol)

Chemical Properties


CAS No. :

Formula: C₂₀H₃₄O₄

Molecular Weight: 338.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Podoblastin A (D-mannitol) is a useful organic compound for research related to life sciences and the catalog number is TSP-41077145.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9543 mL	14.7715 mL	29.543 mL
5 mM	0.5909 mL	2.9543 mL	5.9086 mL
10 mM	0.2954 mL	1.4771 mL	2.9543 mL
50 mM	0.0591 mL	0.2954 mL	0.5909 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481