

## Compound T127967

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>38</sub>O<sub>8</sub>

Molecular Weight: 418.53

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(2S,3R,4S,5S,6R)-2-({2-[(2R,4aR,5R,8R)-5,8-dihydroxy-4a,8-dimethyl-decahydronaphthalen-2-yl]propan-2-yl}oxy)-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3893 mL	11.9466 mL	23.8931 mL
5 mM	0.4779 mL	2.3893 mL	4.7786 mL
10 mM	0.2389 mL	1.1947 mL	2.3893 mL
50 mM	0.0478 mL	0.2389 mL	0.4779 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481