

XIAP degrader-1

Chemical Properties

CAS No. :

Formula: C34H45N5O4

Molecular Weight: 587.75

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	XIAP degrader-1 is a small primary amine molecule that promotes the degradation of X-linked apoptosis inhibitory protein (XIAP).
Targets(IC50)	Others,IAP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7014 mL	8.507 mL	17.014 mL
5 mM	0.3403 mL	1.7014 mL	3.4028 mL
10 mM	0.1701 mL	0.8507 mL	1.7014 mL
50 mM	0.034 mL	0.1701 mL	0.3403 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481