


Schisandrin C epoxide

Chemical Properties

CAS No. :	81345-36-0
Formula:	C ₂₂ H ₂₂ O ₇
Molecular Weight:	398.411
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Schisandrin C epoxide is a natural product for research related to life sciences. The catalog number is TN4958 and the CAS number is 81345-36-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.510 mL	12.5499 mL	25.0998 mL
5 mM	0.502 mL	2.510 mL	5.020 mL
10 mM	0.251 mL	1.255 mL	2.510 mL
50 mM	0.0502 mL	0.251 mL	0.502 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481