


Enalapril-D5 maleate

Chemical Properties

CAS No. :	349554-02-5
Formula:	C ₂₄ H ₃₂ N ₂ O ₉
Molecular Weight:	497.55
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Enalapril-D5 maleate is deuterium labeled Enalapril maleate (TMSM-1101), which is an angiotensin converting enzyme (ACE) inhibitor.
Targets(IC50)	Angiotensin-converting Enzyme (ACE)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0098 mL	10.0492 mL	20.0985 mL
5 mM	0.402 mL	2.0098 mL	4.0197 mL
10 mM	0.201 mL	1.0049 mL	2.0098 mL
50 mM	0.0402 mL	0.201 mL	0.402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481