


ORN-0829 hydrate

Chemical Properties

CAS No. : 2265899-50-9
 Formula: C92H90F4N28O9
 Molecular Weight: 1807.91
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	ORN-0829 hydrate is a potent dual orexin 1/2 receptor antagonist for the treatment of insomnia.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5531 mL	2.7656 mL	5.5312 mL
5 mM	0.1106 mL	0.5531 mL	1.1062 mL
10 mM	0.0553 mL	0.2766 mL	0.5531 mL
50 mM	0.0111 mL	0.0553 mL	0.1106 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481