


## Erigoster B

## Chemical Properties

CAS No. : 849777-61-3  
 Formula: C<sub>26</sub>H<sub>24</sub>O<sub>13</sub>  
 Molecular Weight: 544.46  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Erigoster B is a natural product
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8367 mL	9.1834 mL	18.3668 mL
5 mM	0.3673 mL	1.8367 mL	3.6734 mL
10 mM	0.1837 mL	0.9183 mL	1.8367 mL
50 mM	0.0367 mL	0.1837 mL	0.3673 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481