


Cyanoside F

Chemical Properties

CAS No. :	1800029-50-8
Formula:	C ₄₁ H ₆₂ O ₁₅
Molecular Weight:	794.92
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cyanoside F is a natural product for research related to life sciences. The catalog number is TN5316 and the CAS number is 1800029-50-8.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.258 mL	6.2899 mL	12.5799 mL
5 mM	0.2516 mL	1.258 mL	2.516 mL
10 mM	0.1258 mL	0.629 mL	1.258 mL
50 mM	0.0252 mL	0.1258 mL	0.2516 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481