


2,4-Dihydroxyphenylacetyl-L-asparagine

Chemical Properties

CAS No. :	111872-98-1
Formula:	C ₁₂ H ₁₄ N ₂ O ₆
Molecular Weight:	282.252
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,4-Dihydroxyphenylacetyl-L-asparagine is a natural product for research related to life sciences. The catalog number is T37949 and the CAS number is 111872-98-1.
Targets(IC50)	Others,Antibacterial

Solubility Information

Solubility	H ₂ O: 50 mM,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.543 mL	17.7148 mL	35.4296 mL
5 mM	0.7086 mL	3.543 mL	7.0859 mL
10 mM	0.3543 mL	1.7715 mL	3.543 mL
50 mM	0.0709 mL	0.3543 mL	0.7086 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481