


Xenyhexenic acid

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 964-82-9 |
| Formula: | C ₁₈ H ₁₈ O ₂ |
| Molecular Weight: | 266.34 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|------------------------------------|
| Description | Xenyhexenic acid is a biochemical. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 3.7546 mL | 18.773 mL | 37.546 mL |
| 5 mM | 0.7509 mL | 3.7546 mL | 7.5092 mL |
| 10 mM | 0.3755 mL | 1.8773 mL | 3.7546 mL |
| 50 mM | 0.0751 mL | 0.3755 mL | 0.7509 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481