


## Piperanine

## Chemical Properties

CAS No. :	23512-46-1
Formula:	C17H21NO3
Molecular Weight:	287.359
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Piperanine is a useful organic compound for research related to life sciences. The catalog number is T125252 and the CAS number is 23512-46-1.
Targets(IC50)	Others,Antifungal

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.480 mL	17.3998 mL	34.7996 mL
5 mM	0.696 mL	3.480 mL	6.9599 mL
10 mM	0.348 mL	1.740 mL	3.480 mL
50 mM	0.0696 mL	0.348 mL	0.696 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481