


cis-Miyabenol C

Chemical Properties

CAS No. :	168037-22-7
Formula:	C42H32O9
Molecular Weight:	680.709
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	cis-Miyabenol C exhibits antioxidative activities and shows strong hyaluronidase inhibitory effects. cis-Miyabenol C also has ecdysteroid antagonist activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4691 mL	7.3453 mL	14.6905 mL
5 mM	0.2938 mL	1.4691 mL	2.9381 mL
10 mM	0.1469 mL	0.7345 mL	1.4691 mL
50 mM	0.0294 mL	0.1469 mL	0.2938 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481