

## Sibiriquinone B

## Chemical Properties

CAS No. :	723300-09-2
Formula:	C19H22O2
Molecular Weight:	282.38
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Sibiriquinone B is a useful organic compound for research related to life sciences and the catalog number is T124328.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5413 mL	17.7066 mL	35.4133 mL
5 mM	0.7083 mL	3.5413 mL	7.0827 mL
10 mM	0.3541 mL	1.7707 mL	3.5413 mL
50 mM	0.0708 mL	0.3541 mL	0.7083 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481