


Sirt2-IN-6

Chemical Properties

CAS No. :	2243151-81-5
Formula:	C ₂₆ H ₂₆ N ₆ O ₃ S
Molecular Weight:	502.59
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sirt2-IN-6 is a potent and selective inhibitor of SIRT2 (IC ₅₀ : 0.815 μM) and can be used to study cancer.
Targets(IC ₅₀)	Others,Sirtuin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9897 mL	9.9485 mL	19.8969 mL
5 mM	0.3979 mL	1.9897 mL	3.9794 mL
10 mM	0.199 mL	0.9948 mL	1.9897 mL
50 mM	0.0398 mL	0.199 mL	0.3979 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481