

Compound 0449-0045

## Chemical Properties


CAS No. :

Formula: C42H64O19

Molecular Weight: 872.955

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0449-0045 is a useful organic compound for research related to life sciences and the catalog number is T131673.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1455 mL	5.7276 mL	11.4553 mL
5 mM	0.2291 mL	1.1455 mL	2.2911 mL
10 mM	0.1146 mL	0.5728 mL	1.1455 mL
50 mM	0.0229 mL	0.1146 mL	0.2291 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481