

Phospho-L-Arginine (lithium salt hydrate)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature,Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Phospho-L-Arginine (lithium salt hydrate) is a L-arginine derivative suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Amino Acids and Derivatives

Solubility Information

Solubility	H2O: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481