

N-terminally acetylated Endomorphin-1 acetate

Chemical Properties

CAS No. :


Formula: C38H44N6O8

Molecular Weight: 712.79

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N-terminally acetylated Endomorphin-1 acetate (Ac-L-Tyr-L-Pro-L-Trp-L-Phe-CONH ₂) is a modified Endomorphin-1.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4029 mL	7.0147 mL	14.0294 mL
5 mM	0.2806 mL	1.4029 mL	2.8059 mL
10 mM	0.1403 mL	0.7015 mL	1.4029 mL
50 mM	0.0281 mL	0.1403 mL	0.2806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481