


## 3''-Galloylquercitrin

## Chemical Properties

CAS No. :	503446-90-0
Formula:	C <sub>28</sub> H <sub>24</sub> O <sub>15</sub>
Molecular Weight:	600.48
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1. 3''-Galloylquercitrin possesses the activity for PTK inhibition.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6653 mL	8.3267 mL	16.6533 mL
5 mM	0.3331 mL	1.6653 mL	3.3307 mL
10 mM	0.1665 mL	0.8327 mL	1.6653 mL
50 mM	0.0333 mL	0.1665 mL	0.3331 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481