


5-Iodoarabinouridine

Chemical Properties

CAS No. :	3052-06-0
Formula:	C ₉ H ₁₁ N ₂ O ₆
Molecular Weight:	370.1
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Iodoarabinouridine is a Halo-nucleoside; Arabino-nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.702 mL	13.5099 mL	27.0197 mL
5 mM	0.5404 mL	2.702 mL	5.4039 mL
10 mM	0.2702 mL	1.351 mL	2.702 mL
50 mM	0.054 mL	0.2702 mL	0.5404 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481