


## 2',3'-Bis-O-(t-butyltrimethylsilyl)-5'-O-(4,4'-dimethyltriphenylmethyl)uridine

### Chemical Properties

CAS No. :	82444-76-6
Formula:	C42H58N2O8Si2
Molecular Weight:	775.09
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



### Biological Description

Description	2',3'-Bis-O-(t-butyltrimethylsilyl)-5'-O-(4,4'-dimethyltriphenylmethyl)uridine is a nucleoside Derivative - Other modified nucleoside.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2902 mL	6.4509 mL	12.9017 mL
5 mM	0.258 mL	1.2902 mL	2.5803 mL
10 mM	0.129 mL	0.6451 mL	1.2902 mL
50 mM	0.0258 mL	0.129 mL	0.258 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481