


## Scopadulciol

## Chemical Properties

CAS No. :	136565-26-9
Formula:	C <sub>27</sub> H <sub>36</sub> O <sub>4</sub>
Molecular Weight:	424.57
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Scopadulciol is an enzyme and proton pump inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3553 mL	11.7766 mL	23.5532 mL
5 mM	0.4711 mL	2.3553 mL	4.7106 mL
10 mM	0.2355 mL	1.1777 mL	2.3553 mL
50 mM	0.0471 mL	0.2355 mL	0.4711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481